



Product specification

PCR single tube, 0.5 ml, PCR Performance Tested, colour mix, PP, flat cap



Product information

Product number	72.735.992
----------------	------------

Product descriptions

Name	PCR single tube, 0.5 ml, PCR Performance Tested, colour mix, PP, flat cap
Product number	72.735.992
Product description	PCR single tube, 0.5 ml, PCR Performance Tested, red/green/yellow/blue/violet/orange, material: PP, with anti-contamination shield, flat cap, per colour 500 piece(s)/minigrip bag

Material & colours

Color of product	Colour mix
Profile	High Profile

Size

Diameter	7.9 mm
Width of product	10 mm
Length of product	30.9 mm
Height of product	29.5 mm
Max. volume	500 µl

Sample volume	0.5 ml
---------------	--------

Product characteristics

Base shape	conical
Graduation	yes
Cap	flat cap
Closure type	Push cap
Autoclavable (max. °C)	121 °C
Suitable for qPCR	no
Centrifugation max (RCF)	8000 x g
Centrifugation temperature	22 °C
Centrifugation time	5 Min

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Piece(s) / inner box	3000
Piece(s) / case	3000
Piece(s) / pallet	198000
Height of case	155 mm
Width of case	323 mm
Depth of case	428 mm
Volume of case	0.0214 cbm
Weight of case	1.671 kg
EAN inner box	4038917120067
EAN case	4038917017053

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

