



Product specification

PCR single tube, 0.5 ml, PCR Performance Tested, colour mix, PP, flat cap



Product information

| | |
|----------------|------------|
| Product number | 72.735.992 |
|----------------|------------|

Product descriptions

| | |
|---------------------|--|
| Name | PCR single tube, 0.5 ml, PCR Performance Tested, colour mix, PP, flat cap |
| Product number | 72.735.992 |
| Product description | PCR single tube, 0.5 ml, PCR Performance Tested, red/green/yellow/blue/violet/orange, material: PP, with anti-contamination shield, flat cap, per colour 500 piece(s)/minigrip bag |

Material & colours

| | |
|------------------|--------------------|
| Product material | Polypropylene (PP) |
| Color of product | Colour mix |
| Profile | High Profile |

Size

| | |
|-------------------|---------|
| Diameter | 7.9 mm |
| Width of product | 10 mm |
| Length of product | 30.9 mm |
| Height of product | 29.5 mm |

| | |
|---------------|--------|
| Max. volume | 500 µl |
| Sample volume | 0.5 ml |

Product characteristics

| | |
|----------------------------|----------|
| Base shape | conical |
| Graduation | yes |
| Cap | flat cap |
| Closure type | Push cap |
| Autoclavable (max. °C) | 121 °C |
| Suitable for qPCR | no |
| Centrifugation max (RCF) | 8000 x g |
| Centrifugation temperature | 22 °C |
| Centrifugation time | 5 Min |

Purity & certification

| | |
|------------------|------------------------|
| Product category | diagnostic device, CE |
| CE certified | CE |
| Purity standard | PCR Performance Tested |
| Batched | yes |

Packaging

| | |
|------------------------|---------------|
| Minimum order quantity | 3000 |
| Piece(s) / inner box | 3000 |
| Piece(s) / case | 3000 |
| Piece(s) / pallet | 198000 |
| Height of case | 155 mm |
| Width of case | 323 mm |
| Depth of case | 428 mm |
| Volume of case | 0.0214 cbm |
| Weight of case | 1.671 kg |
| EAN inner box | 4038917120067 |
| EAN case | 4038917017053 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support
and is valid without signature.



SARSTEDT