


Product specification

PCR plate without skirt, 96 well, transparent, High-Profile, 200 µl, Biosphere® plus, PP



Product information

Order number	72.985
--------------	--------

Product descriptions

Name	PCR plate without skirt, 96 well, transparent, High-Profile, 200 µl, Biosphere® plus, PP
Product number	72.985
Product description	PCR plate Multiply®, 96 wells, transp., 200 µl max. vol., Biosphere® plus (sterile, free of DNA, DNase/RNase, PCR inhib., ATP and pyrogens/endotoxins), 12 str. of 8 tubes (72.985.002) preloaded in PCR work tray (95.987.002), 1 pc./paper/plastic pack.

Product characteristics

Number of wells	96
Construction of the edges of the plate	without skirt
Base shape	conical
Cap	flat lid
Type of cap	Push cap
Autoclavable (max. °C)	121 °C
Suitable for qPCR	yes
Centrifugation max (RCF)	3000 x g
Centrifugation temperature	22 °C
Centrifugation time	5 Min

Size

Diameter	5.4 mm
Width of product	79 mm
Length of product	115 mm
Height of product	20.7 mm
Max. volume	200 µl
Sample volume	0.2 ml

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Profile	High-Profile

Purity & certification

Accordance with MDD or IVDD	IVD product
CE certified	CE
Purity standard	Biosphere plus
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Lot or serial number	yes

Packaging

Piece(s) / inner box	20
Piece(s) / case	20
Piece(s) / pallet	3360
Height of case	150 mm
Width of case	218 mm
Depth of case	288 mm
Volume of case	0.0094 cbm
Weight of case	0.76 kg
EAN inner box	4038917084840
EAN case	4038917003612

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.

