



Product specification

PCR single tube, 0.2 ml, PCR Performance Tested, colour mix, PP, flat cap



Product information

Product number	72.737.992
----------------	------------

Product descriptions

Name	PCR single tube, 0.2 ml, PCR Performance Tested, colour mix, PP, flat cap
------	---

Product number	72.737.992
----------------	------------

Product description	PCR single tube, 0.2 ml, PCR Performance Tested, red/green/yellow/blue/violet/orange, material: PP, with anti-contamination shield, flat cap, suitable for qPCR, per colour 500 piece(s)/minigrip bag
---------------------	---

Material & colours

Color of product	Colour mix
------------------	------------

Profile	High Profile
---------	--------------

Size

Diameter	7.5 mm
----------	--------

Width of product	7.5 mm
------------------	--------

Length of product	21.7 mm
-------------------	---------

Height of product	20.65 mm
-------------------	----------

Max. volume	200 µl
-------------	--------

Sample volume	0.2 ml
---------------	--------

Product characteristics

Base shape	conical
Graduation	no
Cap	flat cap
Closure type	Push cap
Autoclavable (max. °C)	121 °C
Suitable for qPCR	yes
Centrifugation max (RCF)	8000 x g
Centrifugation temperature	22 °C
Centrifugation time	5 Min

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Piece(s) / inner box	3000
Piece(s) / case	3000
Piece(s) / pallet	432000
Height of case	145 mm
Width of case	213 mm
Depth of case	323 mm
Volume of case	0.01 cbm
Weight of case	0.759 kg
EAN inner box	4038917120197
EAN case	4038917067676

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

