




Product specification

PCR strip of lids, PCR Performance Tested, PP, highly transparent, flat cap



Product information

Product number	65.989.002
----------------	------------

Product descriptions

Name	PCR strip of lids, PCR Performance Tested, PP, highly transparent, flat cap
Product number	65.989.002
Product description	PCR strip of lids, for PCR plates (REF 72.1978/.010, 72.1979/.010/.003/.201/.700, 72.1980/.010/.201/.600/.700, 72.1981/.010) & PCR strips (REF 72.985.002/.092/.992), PCR Performance Tested, material: PP, highly transparent, flat cap, suitable for qPCR, 120 piece(s)/minigrip bag

Material & colours

Product material	Polypropylene (PP)
Color of product	highly transparent

Size

Diameter	4.47 mm
Width of product	7.3 mm
Length of product	77.2 mm

Height of product	3.25 mm
-------------------	---------

Product characteristics

Cap	flat cap
Closure type	Push cap
Autoclavable (max. °C)	121 °C
Pierceable	no
Suitable for qPCR	yes

Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Purity standard	PCR Performance Tested
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order quantity	480
Piece(s) / inner box	120
Piece(s) / case	480
Piece(s) / pallet	155520
Height of case	100 mm
Width of case	191 mm
Depth of case	268 mm
Volume of case	0.0051 cbm
Weight of case	0.3802 kg
EAN inner box	4038917259910
EAN case	4038917259903

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT