

Product specification

Screw cap, yellow, suitable for protective container 126 x 30 mm



Product information

Product number	65.679.004
----------------	------------

Product descriptions

Name	Screw cap, yellow, suitable for protective container 126 x 30 mm
Product number	65.679.004
Product description	Screw cap, yellow, suitable for protective container 126 x 30 mm, 50 piece(s)/case

Product characteristics

Suitable for	Protective container 126 x 30 mm
--------------	----------------------------------

Size

Length of product	20 mm
-------------------	-------

Material & colours

Product material	Polypropylene (PP)
Color of product	yellow

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Piece(s) / inner box	50
Piece(s) / case	250
Piece(s) / pallet	41250
Height of case	120 mm

Width of case	228 mm
Depth of case	338 mm
Volume of case	0.0092 cbm
Weight of case	1.1 kg
EAN inner box	4038917115834
EAN case	4038917014939

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.

