



Product specification

Scintillation vials



Product information

Product number	73.662.500
----------------	------------

Product descriptions

Name	Scintillation vials
Product number	73.662.500
Product description	Scintillation vials, material: HD-PE, 20 ml, (LxØ): 58 x 27 mm, 500 piece(s)/bag

Material & colours

Closure material	High Density Polyethylene (HD-PE)
------------------	-----------------------------------

Size

Diameter	27 mm
Length of product	58 mm
Height of product	58 mm
Diameter of opening	17 mm
Max. volume	20 ml

Product characteristics

Base shape	Flat base
Graduation	yes
Closure type	Screw cap

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Piece(s) / inner box	500
Piece(s) / case	500
Piece(s) / pallet	21000
Height of inner box	235 mm
Width of inner box	370 mm
Depth of inner box	380 mm
Height of case	235 mm
Width of case	370 mm
Depth of case	380 mm
Volume of case	0.033 cbm
Weight of case	4 kg
EAN inner box	4038917115728
EAN case	4038917014861

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT