



Product specification

Sample tube, suitable for Coulter Counter



Product information

Product number	73.650.500
----------------	------------

Product descriptions

Name	Sample tube, suitable for Coulter Counter
Product number	73.650.500
Product description	Sample tube, suitable for Coulter Counter, cap enclosed, material: PS, (LxØ): 60 x 34 mm, 250 piece(s)/bag

Material & colours

Closure material	Low Density Polyethylene (LD-PE)
------------------	----------------------------------

Size

Diameter	34 mm
Length of product	60 mm
Diameter of opening	12 mm

Product characteristics

Base shape	Flat base
Graduation	no
Closure type	Push cap

Purity & certification

Product category	In-vitro diagnostic (IVD)
------------------	---------------------------

CE certified

CE - manufacturer's self-declaration

Batched

yes

Packaging

Piece(s) / inner box

250

Piece(s) / case

500

Piece(s) / pallet

18000

Height of case

310 mm

Width of case

398 mm

Depth of case

398 mm

Volume of case

0.0491 cbm

Weight of case

4.51 kg

EAN inner box

4038917157841

EAN case

4038917157834

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT