



Product specification

Sample tube, suitable for Coulter Counter



Product information

Product number	73.650.500
----------------	------------

Product descriptions

Name	Sample tube, suitable for Coulter Counter
Product number	73.650.500
Product description	Sample tube, suitable for Coulter Counter, cap enclosed, material: PS, (LxØ): 60 x 34 mm, 250 piece(s)/bag

Material & colours

Product material	Polystyrene (PS)
Closure material	Low Density Polyethylene (LD-PE)

Size

Diameter	34 mm
Length of product	60 mm
Diameter of opening	12 mm

Product characteristics

Base shape	Flat base
Graduation	no
Closure type	Push cap

Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order quantity	500
Piece(s) / inner box	250
Piece(s) / case	500
Piece(s) / pallet	18000
Height of case	300 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0475 cbm
Weight of case	4.51 kg
EAN inner box	4038917157841
EAN case	4038917157834

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT