



Product specification

RIA tube, suitable for gamma counters, transparent



Product information

Product number	73.1055
----------------	---------

Product descriptions

Name	RIA tube, suitable for gamma counters, transparent
Product number	73.1055
Product description	RIA tube, suitable for gamma counters, material: PP, 0.5 ml, (LxØ): 22 x 8 mm, transparent, 1,000 piece(s)/bag

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	8 mm
Length of product	22 mm
Height of product	22 mm
Diameter of opening	29 mm
Max. volume	0.5 ml

Product characteristics

Base shape	Conical base
Graduation	no
Closure type	Push cap

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Piece(s) / inner box	1000
Piece(s) / case	10000
Piece(s) / pallet	480000
Height of case	205 mm
Width of case	378 mm
Depth of case	388 mm
Volume of case	0.0301 cbm
Weight of case	5.05 kg
EAN inner box	4038917085175
EAN case	4038917003865

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT