



## Product specification

### Sample tube, suitable for auto-analyzer, transparent



## Product information

Product number	73.646
----------------	--------

## Product descriptions

Name	Sample tube, suitable for Auto-analyzer, transparent
Product number	73.646
Product description	Sample tube, suitable for auto-analyzer, material: PS, 3.5 ml, (LxØ): 38 x 14 mm, transparent, 1,000 piece(s)/bag

## Material & colours

Product material	Polystyrene (PS)
Color of product	transparent

## Size

Diameter	14 mm
Length of product	38 mm
Height of product	38 mm
Diameter of opening	8 mm
Max. volume	3.5 ml

## Product characteristics

Base shape	Conical base
Graduation	no

Closure type	Push cap
--------------	----------

## Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Batched	yes

## Packaging

Minimum order quantity	5000
Piece(s) / inner box	1000
Piece(s) / case	5000
Piece(s) / pallet	120000
Height of case	415 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0657 cbm
Weight of case	9.62 kg
EAN inner box	4038917083096
EAN case	4038917002066

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT