



Product specification

Sample tube, suitable for auto-analyzer, transparent



Product information

Product number	73.641
----------------	--------

Product descriptions

Name	Sample tube, suitable for Auto-analyzer, transparent
Product number	73.641
Product description	Sample tube, suitable for auto-analyzer, material: PS, 1.5 ml, (LxØ): 24 x 14 mm, transparent, 1,000 piece(s)/bag

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	14 mm
Length of product	24 mm
Height of product	24 mm
Diameter of opening	8 mm
Max. volume	1.5 ml

Product characteristics

Base shape	Conical base
Graduation	no
Closure type	Push cap

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Piece(s) / inner box	1000
Piece(s) / case	5000
Piece(s) / pallet	210000
Height of case	250 mm
Width of case	378 mm
Depth of case	375 mm
Volume of case	0.0354 cbm
Weight of case	5.37 kg
EAN inner box	4038917083072
EAN case	4038917002042

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT