



Product specification

Lid for sample tube 73.666, length: 8 mm, LD-PE



Product information

Product number	65.718
----------------	--------

Product descriptions

Name	Lid for sample tube 73.666, Length: 8 mm, LD-PE
Product number	65.718
Product description	Lid for sample tube 73.666, natural, length: 8 mm, material: LD-PE, 1,000 piece(s)/bag

Material & colours

Product material	Low Density Polyethylene (LD-PE)
Color of product	natural

Size

Length of product	8 mm
-------------------	------

Product characteristics

Cap	Skirted push cap
Closure type	Push cap

Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order quantity	10000
Piece(s) / inner box	1000
Piece(s) / case	10000
Piece(s) / pallet	360000
Height of case	300 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0475 cbm
Weight of case	6.97 kg
EAN inner box	4038917083669
EAN case	4038917002561

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT