



Product specification

Lid for sample tube 73.641 & 73.646, length: 8 mm, LD-PE



Product information

Product number	65.649
----------------	--------

Product descriptions

Name	Lid for sample tube 73.641 & 73.646, length: 8 mm, LD-PE
Product number	65.649
Product description	Lid for sample tube 73.641 & 73.646, natural, length: 8 mm, material: LD-PE, 1,000 piece(s)/bag

Material & colours

Color of product	natural
------------------	---------

Size

Length of product	8 mm
-------------------	------

Product characteristics

Cap	Skirted push cap
-----	------------------

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Piece(s) / inner box	1000
Piece(s) / case	10000
Piece(s) / pallet	420000
Height of inner box	45 mm
Width of inner box	300 mm
Depth of inner box	350 mm
Height of case	250 mm
Width of case	378 mm
Depth of case	375 mm
Volume of case	0.0354 cbm
Weight of case	5.53 kg
EAN inner box	4038917083119
EAN case	4038917002073

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT