



Product specification

Cryobox, 132 x 132 x 53 mm, format: 10 x 10, for 100 collection tubes



Product information

Product number	93.874.610
----------------	------------

Product descriptions

Name	Cryobox, 132 x 132 x 53 mm, format: 10 x 10, for 100 Collection tubes
------	---

Product number	93.874.610
----------------	------------

Product description	Cryobox, with numerical coding at each aperture, for low-temperature storage, material: PC, blue, slip-on lid with ventilation function, cap: transparent, (LxWxH): 132 x 132 x 53 mm, format: 10 x 10, for 100 collection tubes, for CryoPure tubes 1.2–2.0 ml internal thread, 5 piece(s)/bag
---------------------	---

Material & colours

Color of product	blue
Color of cap	transparent

Size

Width of product	132 mm
Length of product	132 mm
Height of product	53 mm

Product characteristics

Design	with numerical coding at each aperture
Graduation	no
Storage capacity	100
Cap	slip-on lid with ventilation function
Closure type	Slip-on lid

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	yes

Packaging

Piece(s) / case	20
Piece(s) / pallet	720
Height of inner box	130 mm
Width of inner box	130 mm
Depth of inner box	267 mm
Height of case	270 mm
Width of case	135 mm
Depth of case	540 mm
Volume of case	0.0197 cbm
Weight of case	3.71 kg
EAN inner box	4038917122931
EAN case	4038917415101

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT