



Product specification

CryoPure tubes, 1.2 ml, QuickSeal screw cap, colour mix



Product information

Product number	72.377.992
----------------	------------

Product descriptions

Name	CryoPure tubes, 1.2 ml, Quickseal screw cap, Colour mix
Product number	72.377.992
Product description	CryoPure tubes, 1.2 ml, tube: PP, QuickSeal screw cap, cap assembled, HD-PE, colour mix, external thread, nominal volume: 1 ml, Cryo Performance Tested, 50 piece(s)/bag

Material & colours

Tube material	Polypropylene (PP)
Closure material	High Density Polyethylene (HD-PE)
Color of cap	Colour mix

Size

Diameter	12.65 mm
Length excluding cap	29.8 mm
Length of product	36.7 mm
Max. volume	1.2 ml
Sample volume	1 ml

Product characteristics

Label/Print	With print
Color of print/label	white
Vessel type	Tube
Thread	External thread
Graduation	yes
Cap	Quickseal screw cap
Closure type	Screw cap

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	Cryo Performance Tested
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / inner box	500
Piece(s) / case	2000
Piece(s) / pallet	104000
Height of inner box	118 mm
Width of inner box	183 mm
Depth of inner box	293 mm
Height of case	135 mm
Width of case	378 mm
Depth of case	598 mm
Volume of case	0.0305 cbm
Weight of case	4.25 kg
EAN inner box	4038917228428
EAN case	4038917228411

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

