



Product specification

Colour-coded inserts, for CryoPure tubes, colour mix



Product information

Product number	65.386.992
----------------	------------

Product descriptions

Name	Colour-coded inserts, for CryoPure tubes, Colour mix
Product number	65.386.992
Product description	Colour-coded inserts, for CryoPure tubes, colour mix, 100 piece(s)/bag

Material & colours

Product material	Polystyrene (PS)
Color of product	Colour mix

Product characteristics

Graduation	no
------------	----

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order quantity	2500
Piece(s) / inner box	100
Piece(s) / case	2500

Piece(s) / pallet	810000
Height of case	100 mm
Width of case	191 mm
Depth of case	268 mm
Volume of case	0.0051 cbm
Weight of case	0.6 kg
EAN inner box	4038917228541
EAN case	4038917228534

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT