



Product specification

Screw cap, natural, suitable for screw cap micro tube

72.733.201



Product information

Product number	65.710
----------------	--------

Product descriptions

Name	Screw cap, natural, suitable for screw cap micro tube 72.733.201
Product number	65.710
Product description	Screw cap, natural, suitable for screw cap micro tube 72.733.201, for volume limitation to 100 µl, 500 piece(s)/bag

Material & colours

Color of product	natural
------------------	---------

Product characteristics

Cap	flat cap
Closure type	Screw cap

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE
Batched	yes

Packaging

Piece(s) / inner box	500
Piece(s) / case	5000
Piece(s) / pallet	220000
Height of inner box	45 mm
Width of inner box	250 mm
Depth of inner box	285 mm
Height of case	155 mm
Width of case	348 mm
Depth of case	568 mm
Volume of case	0.0306 cbm
Weight of case	3.44 kg
EAN inner box	4038917083607
EAN case	4038917080293

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT