



## Product specification

**Push cap, natural, suitable for micro tube Ø 10.8 mm**



### Product information

Product number	65.697
----------------	--------

### Product descriptions

Name	Push cap, natural, suitable for Micro tube Ø 10.8 mm
Product number	65.697
Product description	Push cap, natural, suitable for micro tube Ø 10.8 mm, flat, 1,000 piece(s)/bag

### Material & colours

Color of product	natural
------------------	---------

### Product characteristics

Cap	Standard cap
Closure type	Push cap

### Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

### Packaging

Piece(s) / inner box	1000
----------------------	------

Piece(s) / case	10000
Piece(s) / pallet	660000
Height of inner box	40 mm
Width of inner box	180 mm
Depth of inner box	210 mm
Height of case	155 mm
Width of case	268 mm
Depth of case	438 mm
Volume of case	0.0182 cbm
Weight of case	3.31 kg
EAN inner box	4038917083492
EAN case	4038917002417

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT