



Product specification

Push cap, natural, suitable for micro tube Ø 10.8 mm



Product information

Product number	65.109
----------------	--------

Product descriptions

Name	Push cap, natural, suitable for Micro tube Ø 10.8 mm
Product number	65.109
Product description	Push cap, natural, suitable for micro tube Ø 10.8 mm, with retaining loop, 1,000 piece(s)/bag

Material & colours

Color of product	natural
------------------	---------

Product characteristics

Cap	Retaining loop cap
Closure type	Push cap

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Piece(s) / inner box	1000
----------------------	------

Piece(s) / case	10000
Piece(s) / pallet	640000
Height of inner box	20 mm
Width of inner box	310 mm
Depth of inner box	365 mm
Height of case	215 mm
Width of case	298 mm
Depth of case	398 mm
Volume of case	0.0255 cbm
Weight of case	3.69 kg
EAN inner box	4038917080873
EAN case	4038917000055

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT