




Product specification

**Archiving cap, light blue, suitable for S-Monovette®,
tubes Ø 13-16 mm**



Product information

Product number	65.647.020
----------------	------------

Product descriptions

Name	Archiving cap, light blue, suitable for S-Monovette®, tubes Ø 13-16 mm
Product number	65.647.020
Product description	Archiving cap, with filter, light blue, suitable for S-Monovette®, tubes Ø 13-16 mm, 1,000 piece(s)/bag

Material & colours

Product material	Polypropylene (PP)
Color of product	light blue
Color of cap	light blue

Product characteristics

Filter	yes
--------	-----

Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order quantity	5000
Piece(s) / inner box	1000
Piece(s) / case	5000
Piece(s) / pallet	210000
Height of case	250 mm
Width of case	378 mm
Depth of case	388 mm
Volume of case	0.0367 cbm
Weight of case	5.37 kg
EAN inner box	4038917250085
EAN case	4038917250078

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT