



## Product specification

### Push cap, natural, suitable for tubes Ø 13 mm



## Product information

Product number	65.818
----------------	--------

## Product descriptions

Name	Push cap, natural, suitable for Tubes Ø 13 mm
Product number	65.818
Product description	Push cap, natural, suitable for tubes Ø 13 mm, flat, 1,000 piece(s)/bag

## Material & colours

Color of product	natural
------------------	---------

## Product characteristics

Cap	flat cap
Closure type	Push cap

## Purity & certification

Product category	no medical device   IVD
Batched	no

## Packaging

Piece(s) / inner box	1000
Piece(s) / case	10000
Piece(s) / pallet	640000

Height of inner box	60 mm
Width of inner box	210 mm
Depth of inner box	230 mm
Height of case	215 mm
Width of case	298 mm
Depth of case	398 mm
Volume of case	0.0255 cbm
Weight of case	5.61 kg
EAN inner box	4038917084291
EAN case	4038917003155

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT