



## Product specification

**Push cap, natural, suitable for tubes Ø 15.7 mm**



### Product information

Product number	65.802
----------------	--------

### Product descriptions

Name	Push cap, natural, suitable for Tubes Ø 15.7 mm
Product number	65.802
Product description	Push cap, natural, suitable for tubes Ø 15.7 mm, flat, 500 piece(s)/bag

### Material & colours

Product material	Low Density Polyethylene (LD-PE)
Color of product	natural

### Product characteristics

Cap	flat cap
Closure type	Push cap

### Purity & certification

Product category	no medical device   IVD
Batched	no

### Packaging

Minimum order quantity	10000
Piece(s) / inner box	500

Piece(s) / case	10000
Piece(s) / pallet	360000
Height of case	300 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0475 cbm
Weight of case	6.51 kg
EAN inner box	4038917084185
EAN case	4038917003049

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT