



Product specification

Push cap, natural, suitable for tubes Ø 24.8 mm



Product information

Product number	65.791
----------------	--------

Product descriptions

Name	Push cap, natural, suitable for Tubes Ø 24.8 mm
Product number	65.791
Product description	Push cap, natural, suitable for tubes Ø 24.8 mm, 1,000 piece(s)/bag

Material & colours

Product material	Low Density Polyethylene (LD-PE)
Color of product	natural

Product characteristics

Cap	Standard cap
Closure type	Push cap

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order quantity	2000
Piece(s) / inner box	1000

Piece(s) / case	2000
Piece(s) / pallet	84000
Height of case	235 mm
Width of case	378 mm
Depth of case	388 mm
Volume of case	0.0345 cbm
Weight of case	4.61 kg
EAN inner box	4038917084147
EAN case	4038917002998

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT