



## Product specification

### Screw cap, natural, suitable for tubes 82 x 13 mm



## Product information

Product number	65.726
----------------	--------

## Product descriptions

Name	Screw cap, natural, suitable for tubes 82 x 13 mm
Product number	65.726
Product description	Screw cap, natural, suitable for tubes 82 x 13 mm, 1,000 piece(s)/bag

## Material & colours

Product material	High Density Polyethylene (HD-PE)
Color of product	natural

## Size

Diameter	16.3 mm
Height of product	9 mm

## Product characteristics

Cap	Screw cap with cotton bud holder
Closure type	Screw cap

## Purity & certification

Product category	no medical device   IVD
Latex-free	yes

Batched	yes
---------	-----

## Packaging

Minimum order quantity	5000
Piece(s) / inner box	1000
Piece(s) / case	5000
Piece(s) / pallet	420000
Height of case	125 mm
Width of case	358 mm
Depth of case	388 mm
Volume of case	0.0174 cbm
Weight of case	4.4 kg
EAN inner box	4038917083744
EAN case	4038917002646

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT