



## Product specification

### Ventilation cap, natural, suitable for tubes Ø 16 mm



## Product information

Product number	65.725
----------------	--------

## Product descriptions

Name	Ventilation cap, natural, suitable for Tubes Ø 16 mm
Product number	65.725
Product description	Ventilation cap, natural, suitable for tubes Ø 16 mm, 100 piece(s)/bag

## Material & colours

Color of product	natural
------------------	---------

## Product characteristics

Cap	Ventilation cap
Closure type	Push cap

## Purity & certification

Product category	no medical device   IVD
Batched	no

## Packaging

Piece(s) / inner box	100
Piece(s) / case	5000
Piece(s) / pallet	120000

Height of case	415 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0657 cbm
Weight of case	5.88 kg
EAN inner box	4038917083737
EAN case	4038917002639

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT