



## Product specification

### Ventilation cap, natural, suitable for tubes Ø 12 mm



## Product information

Product number	65.722
----------------	--------

## Product descriptions

Name	Ventilation cap, natural, suitable for Tubes Ø 12 mm
Product number	65.722
Product description	Ventilation cap, natural, suitable for tubes Ø 12 mm, 100 piece(s)/bag

## Material & colours

Product material	Low Density Polyethylene (LD-PE)
Color of product	natural

## Product characteristics

Cap	Ventilation cap
Closure type	Push cap

## Purity & certification

Product category	no medical device   IVD
Latex-free	yes
Batched	yes

## Packaging

Minimum order quantity	5000
------------------------	------

Piece(s) / inner box	100
Piece(s) / case	5000
Piece(s) / pallet	240000
Height of case	215 mm
Width of case	378 mm
Depth of case	388 mm
Volume of case	0.0315 cbm
Weight of case	3.71 kg
EAN inner box	4038917083706
EAN case	4038917002608

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT