



Product specification

Skirted push cap, natural, suitable for tubes Ø 11.5 and 12 mm



Product information

Product number	65.719
----------------	--------

Product descriptions

Name	Skirted push cap, natural, suitable for tubes Ø 11.5 and 12 mm
Product number	65.719
Product description	Skirted push cap, natural, suitable for tubes Ø 11.5 and 12 mm, 1,000 piece(s)/bag

Material & colours

Product material	Low Density Polyethylene (LD-PE)
Color of product	natural

Product characteristics

Cap	Skirted push cap
Closure type	Push cap

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order quantity	5000
Piece(s) / inner box	1000
Piece(s) / case	5000
Piece(s) / pallet	320000
Height of case	215 mm
Width of case	298 mm
Depth of case	398 mm
Volume of case	0.0255 cbm
Weight of case	5.04 kg
EAN inner box	4038917083676
EAN case	4038917002578

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT