




## Product specification

### Screw cap, natural, suitable for tubes Ø 16.5 mm



## Product information

Product number	65.714
----------------	--------

## Product descriptions

Name	Screw cap, natural, suitable for Tubes Ø 16.5 mm
Product number	65.714
Product description	Screw cap, natural, suitable for tubes Ø 16.5 mm, 1,000 piece(s)/bag

## Material & colours

Color of product	natural
O-Ring material	Ethylene propylene diene rubber (EPDM)

## Product characteristics

Cap	tall cap
Closure type	Screw cap

## Purity & certification

Product category	no medical device   IVD
Batched	yes

## Packaging

Piece(s) / inner box	1000
----------------------	------

Piece(s) / case	5000
Piece(s) / pallet	180000
Height of case	265 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.042 cbm
Weight of case	7.14 kg
EAN inner box	4038917083638
EAN case	4038917002530

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT