




Product specification

Archiving screw cap, light blue, suitable for tubes Ø 15.3 mm



Product information

Product number	65.645
----------------	--------

Product descriptions

Name	Archiving screw cap, light blue, suitable for tubes Ø 15.3 mm
Product number	65.645
Product description	Archiving screw cap, light blue, suitable for tubes Ø 15.3 mm, 1,000 piece(s)/bag

Material & colours

Color of product	light blue
Tube material	Polystyrene (PS)
Color of cap	light blue

Product characteristics

Cap	flat cap
Closure type	Screw cap

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Piece(s) / inner box	1000
Piece(s) / case	5000
Piece(s) / pallet	240000
Height of inner box	50 mm
Width of inner box	300 mm
Depth of inner box	300 mm
Height of case	215 mm
Width of case	378 mm
Depth of case	388 mm
Volume of case	0.0315 cbm
Weight of case	5.37 kg
EAN inner box	4038917274470
EAN case	4038917274463

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT