




Product specification

Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PP, with print



Product information

Product number	62.554.502
----------------	------------

Product descriptions

Name	Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PP, with print
------	---

Product number	62.554.502
----------------	------------

Product description	Screw cap tube, working volume: 15 ml, (LxØ): 120 x 17 mm, material: PP, conical base, transparent, screw cap, red, cap assembled, with print, label/print: white/blue, with graduations, free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile, 50 piece(s)/bag
---------------------	---

Material & colours

Color of product	transparent
------------------	-------------

Closure material	High Density Polyethylene (HD-PE)
------------------	-----------------------------------

Color of cap	red
--------------	-----

Size

Diameter	17 mm
----------	-------

Length of product	120 mm
-------------------	--------

Base shape	Conical base
Label/Print	With print
Color of print/label	white/blue
Graduation	yes
Cap	standard screw cap
Closure type	Screw cap
Centrifugation max (RCF)	20000 x g

Purity & certification

Satisfies the requirement	IATA
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	20000
Height of case	165 mm
Width of case	398 mm
Depth of case	598 mm
Volume of case	0.0393 cbm
Weight of case	3.74 kg
EAN inner box	4038917113182
EAN case	4038917013932

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.



SARSTEDT