



Product specification

Screw cap tube, 50 ml, (LxØ): 114 x 28 mm, PP, with print



Product information

Product number	62.547.354
----------------	------------

Product descriptions

Name	Screw cap tube, 50 ml, (LxØ): 114 x 28 mm, PP, with print
Product number	62.547.354
Product description	Screw cap tube, working volume: 50 ml, (LxØ): 114 x 28 mm, material: PP, conical base, brown, HD-PE, cap: red, cap assembled, with print, label/print: blue/white, with graduations, free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile, 25 piece(s)/bag

Material & colours

Color of product	brown
Closure material	High Density Polyethylene (HD-PE)
Color of cap	red

Size

Diameter	28 mm
Length of product	114 mm

Product characteristics

Base shape	Conical base
------------	--------------

Label/Print	With print
Color of print/label	blue/white
Graduation	yes
Closure type	Screw cap
Centrifugation max (RCF)	20000 x g

Purity & certification

Satisfies the requirement	ADR
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / inner box	25
Piece(s) / case	300
Piece(s) / pallet	8400
Height of case	255 mm
Width of case	400 mm
Depth of case	590 mm
Volume of case	0.0602 cbm
Weight of case	4.5399 kg
EAN inner box	4038917157759
EAN case	4038917157742

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT