




## Product specification

**Tube, 12 ml, (LxØ): 110 x 17 mm, PS**



## Product information

Product number	57.462.015
----------------	------------

## Product descriptions

Name	Tube, 12 ml, (LxØ): 110 x 17 mm, PS
Product number	57.462.015
Product description	Tube, working volume: 12 ml, (LxØ): 110 x 17 mm, material: PS, conical base, transparent, push cap, with graduations, 250 piece(s)/StackPack

## Material & colours

Color of product	transparent
------------------	-------------

## Size

Diameter	17 mm
Length of product	110 mm

## Product characteristics

Base shape	Conical base
Graduation	yes
Closure type	Push cap
Centrifugation max (RCF)	2500 x g

## Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
---------------------------	--------------------------

Product category

In-vitro diagnostic (IVD)

CE certified

CE - manufacturer's self-declaration

Batched

yes

## Packaging

Piece(s) / inner box

250

Piece(s) / case

1000

Piece(s) / pallet

42000

Height of inner box

201 mm

Width of inner box

114 mm

Depth of inner box

366 mm

Height of case

250 mm

Width of case

378 mm

Depth of case

375 mm

Volume of case

0.0354 cbm

Weight of case

6.96 kg

EAN inner box

4038917381079

EAN case

4038917381062

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT