




## Product specification

**Tube, 12 ml, (LxØ): 110 x 17 mm, PS**



## Product information

Product number	57.462.015
----------------	------------

## Product descriptions

Name	Tube, 12 ml, (LxØ): 110 x 17 mm, PS
Product number	57.462.015
Product description	Tube, working volume: 12 ml, (LxØ): 110 x 17 mm, material: PS, conical base, transparent, push cap, with graduations, 250 piece(s)/StackPack

## Material & colours

Product material	Polystyrene (PS)
Color of product	transparent

## Size

Diameter	17 mm
Length of product	110 mm

## Product characteristics

Base shape	Conical base
Graduation	yes
Closure type	Push cap
Centrifugation max (RCF)	2500 x g

## Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	diagnostic device, CE
CE certified	CE
Batched	yes

## Packaging

Minimum order quantity	1000
Piece(s) / inner box	250
Piece(s) / case	1000
Piece(s) / pallet	42000
Height of inner box	201 mm
Width of inner box	114 mm
Depth of inner box	366 mm
Height of case	250 mm
Width of case	378 mm
Depth of case	388 mm
Volume of case	0.0367 cbm
Weight of case	6.96 kg
EAN inner box	4038917381079
EAN case	4038917381062

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT