



Product specification

Screw cap tube, 50 ml, (LxØ): 115 x 28 mm, PP, with print



Product information

Product number	62.559.205
----------------	------------

Product descriptions

Name	Screw cap tube, 50 ml, (LxØ): 115 x 28 mm, PP, with print
------	---

Product number	62.559.205
----------------	------------

Product description	Screw cap tube, working volume: 50 ml, (LxØ): 115 x 28 mm, material: PP, skirted conical base, transparent, screw cap, red, cap assembled, with print, label/print: white/blue, with graduations, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 25 piece(s)/bag
---------------------	---

Material & colours

Product material	Polypropylene (PP)
------------------	--------------------

Color of product	transparent
------------------	-------------

Closure material	High Density Polyethylene (HD-PE)
------------------	-----------------------------------

Color of cap	red
--------------	-----

Size

Diameter	28 mm
----------	-------

Length of product	115 mm
-------------------	--------

Product characteristics

Base shape	skirted conical base
Label/Print	With print
Color of print/label	white/blue
Graduation	yes
Closure type	Screw cap
Centrifugation max (RCF)	8000 x g

Purity & certification

Satisfies the requirement	IATA
Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order quantity	300
Piece(s) / inner box	25
Piece(s) / case	300
Piece(s) / pallet	8400
Volume of case	0.0649 cbm
Weight of case	5.19 kg
EAN inner box	4038917376402
EAN case	4038917376396

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT