




## Product specification

**Screw cap tube, 4.5 ml, (LxØ): 75 x 12 mm, PP, with print**



## Product information

Product number	62.557.012
----------------	------------

## Product descriptions

Name	Screw cap tube, 4.5 ml, (LxØ): 75 x 12 mm, PP, with print
------	---

Product number	62.557.012
----------------	------------

Product description	Screw cap tube, working volume: 4.5 ml, (LxØ): 75 x 12 mm, material: PP, round base, transparent, screw cap, natural, cap assembled, with print, label/print: white, with graduations, sterile, 100 piece(s)/bag
---------------------	--

## Material & colours

Color of product	transparent
------------------	-------------

Closure material	High Density Polyethylene (HD-PE)
------------------	-----------------------------------

Color of cap	natural
--------------	---------

## Size

Diameter	12 mm
----------	-------

Length of product	75 mm
-------------------	-------

## Product characteristics

Base shape	Round base
------------	------------

Label/Print	With print
-------------	------------

Color of print/label	white
----------------------	-------

Graduation	yes
Closure type	Screw cap
Centrifugation max (RCF)	5000 x g

## Purity & certification

Satisfies the requirement	IATA
Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

## Packaging

Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	56000
Volume of case	0.029 cbm
Weight of case	5.15 kg
EAN inner box	4038917225762
EAN case	4038917225755

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT