



Product specification

Screw cap tube, 4.5 ml, (LxØ): 75 x 12 mm, PP, with print



Product information

| | |
|----------------|------------|
| Product number | 62.557.002 |
|----------------|------------|

Product descriptions

| | |
|---------------------|--|
| Name | Screw cap tube, 4.5 ml, (LxØ): 75 x 12 mm, PP, with print |
| Product number | 62.557.002 |
| Product description | Screw cap tube, working volume: 4.5 ml, (LxØ): 75 x 12 mm, material: PP, round base, transparent, screw cap, natural, cap enclosed, with print, label/print: white, with graduations, 1,000 piece(s)/bag |

Material & colours

| | |
|------------------|-----------------------------------|
| Color of product | transparent |
| Tube material | Polypropylene (PP) |
| Closure material | High Density Polyethylene (HD-PE) |
| Color of cap | natural |

Size

| | |
|-------------------|-------|
| Diameter | 12 mm |
| Length of product | 75 mm |

Product characteristics

| | |
|-------------|------------|
| Base shape | Round base |
| Label/Print | With print |

| | |
|--------------------------|-----------|
| Color of print/label | white |
| Graduation | yes |
| Closure type | Screw cap |
| Autoclavable (max. °C) | 121 °C |
| Centrifugation max (RCF) | 4000 x g |

Purity & certification

| | |
|---------------------------|-------------------------|
| Satisfies the requirement | IATA |
| Product category | no medical device IVD |
| Batched | yes |

Packaging

| | |
|----------------------|---------------|
| Piece(s) / inner box | 1000 |
| Piece(s) / case | 1000 |
| Piece(s) / pallet | 56000 |
| Volume of case | 0.029 cbm |
| Weight of case | 5.15 kg |
| EAN inner box | 4038917113342 |
| EAN case | 4038917065092 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT