



Product specification

Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PP, with print



Product information

Product number	62.554.205
----------------	------------

Product descriptions

Name	Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PP, with print
Product number	62.554.205
Product description	Screw cap tube, working volume: 15 ml, (LxØ): 120 x 17 mm, material: PP, conical base, transparent, screw cap, red, cap assembled, with print, label/print: white/blue, with graduations, free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile, 50 piece(s)/bag

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	red

Size

Diameter	17 mm
Length of product	120 mm

Product characteristics

Base shape	Conical base
Label/Print	With print
Color of print/label	white/blue
Graduation	yes
Cap	standard screw cap
Closure type	Screw cap
Centrifugation max (RCF)	20000 x g

Purity & certification

Satisfies the requirement	IATA
Product category	no medical device IVD
Purity standard	free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order quantity	500
Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	24000
Height of inner box	100 mm
Width of inner box	340 mm
Depth of inner box	300 mm
Height of case	219 mm
Width of case	391 mm
Depth of case	429 mm
Volume of case	0.0366 cbm
Weight of case	4 kg
EAN inner box	4038917113144
EAN case	4038917033497

This document was prepared by EDP support
and is valid without signature.



SARSTEDT