



## Product specification

**Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PP, with print**



## Product information

Product number	62.554.100
----------------	------------

## Product descriptions

Name	Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PP, with print
Product number	62.554.100
Product description	Screw cap tube, working volume: 15 ml, (LxØ): 120 x 17 mm, material: PP, conical base, transparent, screw cap, red, cap assembled, with print, label/print: white/blue, with graduations, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 50 piece(s)/rack

## Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	red

## Size

Diameter	17 mm
Length of product	120 mm

Base shape	Conical base
Label/Print	With print
Color of print/label	white/blue
Graduation	yes
Cap	standard screw cap
Closure type	Screw cap
Centrifugation max (RCF)	20000 x g

## Purity & certification

Satisfies the requirement	IATA
Product category	no medical device   IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

## Packaging

Minimum order quantity	25000
Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	18000
Volume of case	0.0436 cbm
Weight of case	4.56 kg
EAN inner box	4038917376341
EAN case	4038917376334

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT