




Product specification

Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PS, with print



Product information

Product number	62.553.002
----------------	------------

Product descriptions

Name	Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PS, with print
------	---

Product number	62.553.002
----------------	------------

Product description	Screw cap tube, working volume: 15 ml, (LxØ): 120 x 17 mm, material: PS, conical base, transparent, screw cap, red, cap assembled, with print, label/print: white/blue, with graduations, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 50 piece(s)/rack
---------------------	--

Material & colours

Product material	Polystyrene (PS)
Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	red

Size

Diameter	17 mm
Length of product	120 mm

Product characteristics

Base shape	Conical base
Label/Print	With print
Color of print/label	white/blue
Graduation	yes
Closure type	Screw cap
Centrifugation max (RCF)	4100 x g

Purity & certification

Satisfies the requirement	IATA
Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Batched	yes

Packaging

Minimum order quantity	500
Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	18000
Volume of case	0.0436 cbm
Weight of case	5.31 kg
EAN inner box	4038917112949
EAN case	4038917033435

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT