



## Product specification

**Screw cap tube, 7 ml, (LxØ): 82 x 13 mm, PP**



### Product information

Product number	60.550.049
----------------	------------

### Product descriptions

Name	Screw cap tube, 7 ml, (LxØ): 82 x 13 mm, PP
Product number	60.550.049
Product description	Screw cap tube, working volume: 7 ml, (LxØ): 82 x 13 mm, material: PP, round base, transparent, screw cap, natural, cap assembled, sterile, 100 piece(s)/bag

### Material & colours

Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	natural

### Size

Diameter	13 mm
Length of product	82 mm

### Product characteristics

Base shape	Round base
Graduation	no
Closure type	Screw cap

Centrifugation max (RCF)

4000 x g

## Purity & certification

Satisfies the requirement	IATA
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile
Batched	yes

## Packaging

Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	36000
Volume of case	0.0436 cbm
Weight of case	4.1 kg
EAN inner box	4038917112529
EAN case	4038917047340

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT