



Product specification

Screw cap tube, 7 ml, (LxØ): 82 x 13 mm, PP



Product information

Product number	60.550.009
----------------	------------

Product descriptions

Name	Screw cap tube, 7 ml, (LxØ): 82 x 13 mm, PP
Product number	60.550.009
Product description	Screw cap tube, working volume: 7 ml, (LxØ): 82 x 13 mm, material: PP, flat base, transparent, screw cap, natural, cap assembled, sterile, 100 piece(s)/bag

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	natural

Size

Diameter	13 mm
Length of product	82 mm

Product characteristics

Base shape	Flat base
Graduation	no

Closure type	Screw cap
Centrifugation max (RCF)	4000 x g

Purity & certification

Satisfies the requirement	IATA
Product category	no medical device IVD
Purity standard	sterile
Batched	yes

Packaging

Minimum order quantity	50000
Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	36000
Volume of case	0.0436 cbm
Weight of case	4.9 kg
EAN inner box	4038917112246
EAN case	4038917032520

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT