




Product specification

Screw cap tube, 50 ml, (LxØ): 114 x 28 mm, PP, with print



Product information

Product number	62.547.205
----------------	------------

Product descriptions

Name	Screw cap tube, 50 ml, (LxØ): 114 x 28 mm, PP, with print
Product number	62.547.205
Product description	Screw cap tube, working volume: 50 ml, (LxØ): 114 x 28 mm, material: PP, conical base, transparent, HD-PE, cap: red, cap assembled, with print, label/print: white/blue, with graduations, free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile, 25 piece(s)/bag

Material & colours

Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	red

Size

Diameter	28 mm
Length of product	114 mm

Base shape	Conical base
Label/Print	With print
Color of print/label	white/blue
Graduation	yes
Cap	standard screw cap
Closure type	Screw cap
Centrifugation max (RCF)	20000 x g

Purity & certification

Satisfies the requirement	ADR
Product category	no medical device IVD
Purity standard	free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / inner box	25
Piece(s) / case	250
Piece(s) / pallet	9000
Volume of case	0.0492 cbm
Weight of case	4.15 kg
EAN inner box	4038917111881
EAN case	4038917033275

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT