




Product specification

Screw cap tube, 50 ml, (LxØ): 114 x 28 mm, PP, with print



Product information

Product number	62.547.100
----------------	------------

Product descriptions

Name	Screw cap tube, 50 ml, (LxØ): 114 x 28 mm, PP, with print
Product number	62.547.100
Product description	Screw cap tube, working volume: 50 ml, (LxØ): 114 x 28 mm, material: PP, conical base, transparent, HD-PE, cap: red, cap assembled, with print, label/print: white/blue, with graduations, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 25 piece(s)/rack

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	red

Size

Diameter	28 mm
Length of product	114 mm

Product characteristics

Base shape	Conical base
Label/Print	With print
Color of print/label	white/blue
Graduation	yes
Closure type	Screw cap
Centrifugation max (RCF)	20000 x g

Purity & certification

Satisfies the requirement	IATA
Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order quantity	200
Piece(s) / inner box	25
Piece(s) / case	200
Piece(s) / pallet	7200
Volume of case	0.0436 cbm
Weight of case	4.01 kg
EAN inner box	4038917376259
EAN case	4038917376242

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT