



Product specification

Screw cap tube, 8 ml, (LxØ): 57 x 16.5 mm, PP, with print



Product information

Product number	62.542.026
----------------	------------

Product descriptions

Name	Screw cap tube, 8 ml, (LxØ): 57 x 16.5 mm, PP, with print
Product number	62.542.026
Product description	Screw cap tube, working volume: 8 ml, (LxØ): 57 x 16.5 mm, material: PP, flat base, transparent, screw cap, without cap, with print, label/print: white, with graduations, 1,000 piece(s)/bag

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	16.5 mm
Length of product	57 mm
Height of product	57 mm

Product characteristics

Base shape	Flat base
Label/Print	With print
Color of print/label	white
Graduation	yes

Closure type	Screw cap
Autoclavable (max. °C)	121 °C
Centrifugation max (RCF)	4000 x g

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Piece(s) / inner box	1000
Piece(s) / case	1000
Piece(s) / pallet	60000
Volume of case	0.0282 cbm
Weight of case	5.15 kg
EAN inner box	4038917111218
EAN case	4038917033206

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT