

# Product specification

## Screw cap tube, 13 ml, (LxØ): 101 x 16.5 mm, PP



### Product information

Product number	60.541.385
----------------	------------

### Product descriptions

Name	Screw cap tube, 13 ml, (LxØ): 101 x 16.5 mm, PP
Product number	60.541.385
Product description	Screw cap tube, working volume: 13 ml, (LxØ): 101 x 16.5 mm, material: PP, round base with skirted base, transparent, HD-PE, cap: natural, cap enclosed, 500 piece(s)/bag

### Product characteristics

Base shape	Round base with skirted base
Graduation	no
Closure type	Screw cap
Autoclavable (max. °C)	121 °C
Centrifugation max (RCF)	4000 x g

### Size

Diameter	16.5 mm
Length of product	101 mm

### Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	natural

### Purity & certification

Satisfies the requirement	IATA
---------------------------	------

Product category	no medical device   IVD
Batched	yes

## Packaging

Piece(s) / inner box	500
Piece(s) / case	1000
Piece(s) / pallet	28000
Volume of case	0.0649 cbm
Weight of case	7.26 kg
EAN inner box	4038917367790
EAN case	4038917367783

**This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.**

**This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.**

This document was prepared by EDP support and is valid without signature.

