



Product specification

Screw cap tube, 13 ml, (LxØ): 101 x 16.5 mm, PP



Product information

Product number	60.540.352
----------------	------------

Product descriptions

Name	Screw cap tube, 13 ml, (LxØ): 101 x 16.5 mm, PP
Product number	60.540.352
Product description	Screw cap tube, working volume: 13 ml, (LxØ): 101 x 16.5 mm, material: PP, round base, transparent, screw cap, blue, cap assembled, sterile, 500 piece(s)/bag

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	blue

Size

Diameter	16.5 mm
Length of product	101 mm

Product characteristics

Base shape	Round base
Graduation	no

Closure type	Screw cap
Centrifugation max (RCF)	4000 x g

Purity & certification

Satisfies the requirement	IATA
Product category	no medical device IVD
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order quantity	500
Piece(s) / inner box	500
Piece(s) / case	500
Piece(s) / pallet	28000
Volume of case	0.029 cbm
Weight of case	3.45 kg
EAN inner box	4038917367646
EAN case	4038917367639

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT