

Product specification

Screw cap tube, 13 ml, (LxØ): 101 x 16 mm, PP



Product information

Product number	60.540.316
----------------	------------

Product descriptions

Name	Screw cap tube, 13 ml, (LxØ): 101 x 16 mm, PP
Product number	60.540.316
Product description	Screw cap tube, working volume: 13 ml, (LxØ): 101 x 16 mm, material: PP, round base, transparent, PP, cap: natural, cap assembled, sterile, 500 piece(s)/bag

Product characteristics

Base shape	Round base
Graduation	no
Type of cap	Screw cap
Centrifugation max (RCF)	4000 x g

Size

Diameter	16 mm
Length of product	101 mm

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
O-Ring material	Ethylene propylene diene rubber (EPDM)
Cap material	Polypropylene (PP)
Color of cap	natural

Purity & certification

Satisfies the requirement	IATA
---------------------------	------

Product category	no medical device IVD
Purity standard	sterile
Batched	yes

Packaging

Piece(s) / inner box	500
Piece(s) / case	500
Piece(s) / pallet	28000
Volume of case	0.029 cbm
Weight of case	3.43 kg
EAN inner box	4038917375986
EAN case	4038917375979

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.

