



## Product specification

**Screw cap tube, 13 ml, (LxØ): 101 x 16 mm, PP**



### Product information

Product number	60.540.316
----------------	------------

### Product descriptions

Name	Screw cap tube, 13 ml, (LxØ): 101 x 16 mm, PP
Product number	60.540.316
Product description	Screw cap tube, working volume: 13 ml, (LxØ): 101 x 16 mm, material: PP, round base, transparent, screw cap, natural, cap assembled, sterile, 500 piece(s)/bag

### Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
O-Ring material	Ethylene propylene diene rubber (EPDM)
Closure material	Polypropylene (PP)
Color of cap	natural

### Size

Diameter	16 mm
Length of product	101 mm

### Product characteristics

Base shape	Round base
------------	------------

Graduation	no
Closure type	Screw cap
Centrifugation max (RCF)	4000 x g

## Purity & certification

Satisfies the requirement	IATA
Product category	no medical device   IVD
Purity standard	sterile
Batched	yes

## Packaging

Minimum order quantity	500
Piece(s) / inner box	500
Piece(s) / case	500
Piece(s) / pallet	28000
Volume of case	0.029 cbm
Weight of case	3.43 kg
EAN inner box	4038917375986
EAN case	4038917375979

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT