



Product specification

Tube, 7 ml, (LxØ): 50 x 16 mm, PP



Product information

Product number	58.536.300
----------------	------------

Product descriptions

Name	Tube, 7 ml, (LxØ): 50 x 16 mm, PP
Product number	58.536.300
Product description	Tube, working volume: 7 ml, (LxØ): 50 x 16 mm, material: PP, flat base, transparent, push cap, 1,000 piece(s)/bag

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	16 mm
Length of product	50 mm
Height of product	50 mm

Product characteristics

Base shape	Flat base
Graduation	no
Closure type	Push cap
Autoclavable (max. °C)	121 °C

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Piece(s) / inner box	1000
Piece(s) / case	2000
Piece(s) / pallet	72000
Volume of case	0.0492 cbm
Weight of case	5.24 kg
EAN inner box	4038917367554
EAN case	4038917367547

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT